

FTIR instructions (09/09/2016)

Before starting your session:

- 1) If using non-standard detector or source, contact staff at least 1 day before your session

Standard detectors:

Bench: DLATGS and MCT standard
Microscope: MCT and InSb standard


Standard sources:

Glowbar (mid-IR)
Tungsten lamp (Near-IR)

- 2) Check that the instrument is under vacuum (if not under vacuum notify staff and record in logbook)
- 3) Check that the N2 cylinder is not empty (if < 500 psi contact staff)
- 4) Check that the correct source is installed and on
- 5) Check that the correct beamsplitter is setup



You are NEVER allowed to change the beamsplitted



- a. Vent the chamber
 - i. Upper left console (diagram above) use the arrow keys to move between lines and enter to select
 - ii. From "START MENU"
 1. Select VACUUM CONTROL
 - a. SAMPLE CHANGE MENU
 - i. VENT OPTICS
- b. Lift the round lid to check the beamsplitter handle color (**Do not touch the beamsplitter**)

Red handle	KBr	mid-IR
Brown handle	CaF2	Near-IR
Silver	6um Mylar	Far-IR

If you need a beamplitter other than the one installed please contact the staff.

- 6) Check that the correct module is installed.
 - a. Diamond-ATR
 - b. Ge-ATR
 - c. Straight through transmission
 - d. SEAGULL variable angle reflectance
- 7) Evacuate the chamber
 - a. On the console menu:
Select: EVACUATE OPTICS

You are now ready to start your session:

- 1) Start the kiosk session using any internet connected device
 - a. Automatically opens the N2 gas valve to the spectrometer
 - b. You should now see blinking green LED in the front of the instrument
- 2) Start OPUS 4.0 Software:
 - a. If computer not on: log into the computer: (no password – just hit return)
 - b. log into OPUS 4.0 (password is OPUS)
 - c. click return to get to the software
- 3) Click the green test tube icon (advanced Data Collection)
 - a. “Basic” tab:
 - i. Load your experiment file (typically setup during training)
 - ii. Change sample name and other information if needed
 - b. “Advanced” tab
 - i. change filename
 - ii. change scan times if necessary
of scans for background should be at least as long as # of scans for sample
 - iii. change result spectrum if necessary
 - c. “check Signal” Tab **DO NOT CLICK ANYWHERE**
 - i. wait for the signal to appear **DO NOT CLICK ANYWHERE** until the signal appears of the software will likely crash.
 - ii. Once you see signal the intensity (can take a couple of minutes) Compare the expected intensity to your observed intensity:




MIR + DLATGS + KBr + 8mm aperture	23,000
MIR + DLATGS + KBr + 8mm aperture+d-ATR	2,500
MIR + MCT + KBr + 2 mm aperture	26,000
MIR + MCT + KBr + 3 mm aperture + d-ATR	12,600
MIR + DLATGS + CaF2 + 8mm aperture	14,000
NIR + DLATGS + CaF2 + 8mm	16,000
NIR + MCT + CaF2 + 1 mm	19,000
NIR + DLATGS + CaF2 + 8mm + d-ATR	1,400
NIR + MCT + CaF2 + 8 mm + d-ATR	14,000

If using MCT: Make sure you filled up the MCT dewar with LN2

- iii. If intensity is very different check your detector, beamsplitter and aperture in the “optics” tab if correct, contact staff
 - iv.
 - d. “Basic” Tab
 - i. make sure you have your reference in place and the system is under vacuum
Note: with d-ATR reference typically with nothing and crew off the diamond
click “Background Single Channel”
 - ii. Wait for the system to finish collecting the background
 - iii. Vent system
 - iv. Place your sample
 - v. Evacuate System
 - vi. Wait for pressure / signal to stabilize
 - vii. Click “Sample Single Channel”
 - viii. Wait for the system to finish measuring
- 4) You should now see your sample’s the pattern

To export your data to ASCII

- 1) Select the pattern to export typically the first box under the file name (also recommended to export the S and R patterns)
- 2) Click in the menu: File / Save File As
- 3) A window popped up
 - a. "Mode" tab
select: "Data Point Table"
 - b. "Data Point Table" tab
don't check the tab box and write a "," between the two numbers



- c. "Select File" tab
check that the Path is correct and type the file Name for the exported file:
e.g. *MyData.csv*
- 4) repeat for each spectrum

When done:

- 1) **Evacuate the bench**
- 2) **If using a source other than MIR turn it off**
- 3) Close the software
- 4) End your kiosk session (the bench LED will typically revert back to red)