

# Precession Electron Diffraction Detection and Phase Mapping of Retained Austenite and Carbides in a Heat Treated Low Alloy Carbon Steel Using a JEOL ARM 200 TEM with an AppFive Precession Electron Diffraction System

John Mardinly<sup>1</sup>, Amith Darbal<sup>2</sup>

<sup>1</sup> LeRoy Eyring Center for Solid State Science, Arizona State University, Tempe, AZ USA, <sup>2</sup> AppFive LLC, Tempe, AZ USA.

Retained austenite in heat treated steels is usually measured using X-Ray diffraction, with protocols outlined in ASTM Standard E975 [1]. Retained austenite has been measured in high alloy steels with large grains of retained austenite using EBSD [2]. However, low alloy steels heat treated to reduce retained austenite may have small grains of retained austenite that cannot reliably be detected by EBSD [3] or X-Ray diffraction.

A steel alloy of composition: 0.5% Carbon, 1.12% Chromium, 1.6% Tungsten, 1% Silicon, 0.35% Manganese, 0.15% Vanadium and balance Iron was submitted to our laboratory for analysis. The specimen was fully austenitized at 900C, quenched to room temperature then to -200C at which temperature it was held for 48 hours. Finally, it was annealed at 225C for three hours. Retained austenite was measured using X-ray Diffraction using protocols outlined in ASTM Standard E975 [1]. There was no retained austenite detected. A slice of the steel was prepared for TEM by standard dimple and ion milling techniques.

Data acquisition was in a JEOL ARM 200 TEM with a Schottky Field Emission source operated at 200KeV in TEM mode, with an AppFive Precession Electron Diffraction system controlling the beam. The microscope was operated in TEM alpha-L mode, spot size 5 using a 10 micron condenser aperture. The microscope was operated in TEM mode, the CEOS probe corrector did not influence the probe size. This choice of optical modes was driven by the desire to use a highly parallel beam for electron diffraction, rather than the highly convergent beam that results from using the probe corrector in STEM mode. The electron diffraction patterns were collected using a Sting Ray Firewire camera focused onto the fluorescent screen of the microscope. This camera collected patterns 144x144 pixels, in 0.01 seconds per pattern. The field of view was 292x301 nanometers, 92,700 electron diffraction patterns were collected from the field of view in 15 minutes. Each pixel was approximately one nanometer in size.

Post acquisition analysis was completed using Astar data analysis software which uses the template method of Rauch et al [3]. Electron diffraction patterns were calculated for every zone axis for ferrite, austenite and all five carbides listed in ASTM Standard E975 [1]. Matching the templates from 7 phases to the 92,700 patterns acquired is very computationally intensive, and took over 5 hours on a 3.1 GHz Core 7 CPU with 64 bit Windows 7, 16GB of RAM and a 250GB solid state drive.

